Physics Lab-III

Course Code: BEB352 Contact Hours: 15

Credit: 01 (L-0, T-0, P-2)

MM: 25

Course Outline:

Each student has to perform eight experiments, selecting a minimum of three experiments from each of the following group;

- 1. To determine the Frequency of an Electrically Maintained Tuning Fork by Melde_s Experiment.
- 2. To determine the frequency of AC mains by Sonometer.
- 3. To study the response curve of a Series LCR circuit and determine its (a) Resonant frequency, (b) Impedance at resonance and (c) Quality factor Q, and (d) Band width.
- 4. To study the response curve of a Parallel LCR circuit and determine its (a) Anti resonant frequency and (b) Quality factor O.
- 5. To determine the specific resistance of a given wire using Carey Foster_s bridge. 6.

To determine a Low Resistance by a Potentiometer.

wills University

- 7. To calibrate the given ammeter and voltmeter by potentiometer.
- 8. To draw hysteresis curve of a given sample of ferromagnetic material and from this to determine magnetic susceptibility and permeability of the given specimen
- 9. To determine the ballistic constant of a ballistic galvanometer.
- 10. To study the variation of magnetic field along the axis of current carrying circular coil and then to estimate the radius of the coil.
- 11. To determine Self Inductance of a Coil by Anderson_s Bridge using AC

Suggested Reading:

☐ GeetaSanon, BSc Practical Physics, 1stEdn. (2007), R. Chand & Co.

- ☐ B. L. Worsnop and H. T. Flint, Advanced Practical Physics, Asia Publishing House, New Delhi.
- ☐ Indu Prakash and Ramakrishna, A Text Book of Practical Physics Vol 1 & Vol 2, KitabMahal, New Delhi.
- D. P. Khandelwal, A Laboratory Manual of Physics for Undergraduate Classes, Vani Publication House, New Delhi

departues of Education Fascity of Education & Visconium L

layer to Proposity, ProBy 100

Invertis University

Bareilly